

The Chair-Check[®] System

Primary applications: Chair, wheelchair or toilet (Battery powered, portable unit)



*Chair-Check II
Control Unit*

Allow your fall risk individuals the freedom of sitting without personal supervision or restraints, while you maintain constant monitoring with the Chair-Check System: a Chair-Check II Control Unit and a pressure sensitive Sensormat[®] designed for chair use. When the Chair-Check System alarms, an audible signal sounds to alert your staff.

Portability – Uses a common 9V alkaline battery. Battery low charge alert.

Reduce battery changes – A rechargeable battery is available, as is a battery charger/AC adaptor.

Multiple applications – May also monitor a bed or toilet. Functions with any Sensormat.

Easily distinguishable alarm - Four settable alarm tones, one with lower volume.

Simple one button control – Press Reset to put the unit on Hold or to silence an alarm.

Multiple delay settings – 0-5 seconds. Zero, or instant, is best for seated monitoring, because the individual is postured for quick stand up. Longer delay may be set, particularly for bed applications.

Reduced caregiver intervention - The alarm automatically ceases and monitoring resumes instantly when/if

the individual sits back down.

Bathroom privacy – By using a Potty-Check[®] Sensormat, the patient/resident can be monitored while seated on the toilet, allowing the caregiver to stand outside the door. The Potty-Check is adhered to the underside of the seat ring.

Mounting choices – Integral stainless steel clip, or optional wire holster for hanging or wall mounting.

*Potty-Check
Disposable
Sensormat*



*IT Chair Sensormat
(30 day warranty, disposable)*

*LT Chair Sensormat
(90 day warranty, disposable)*



IT Bed Sensormat (30 day warranty, disposable)

LT Bed Sensormat (90 day warranty, disposable)