

Rechargeable Replacement Battery for Chair-Check® II

- ✓ Nickel-metal hydride battery allows recharging over 1,000 times.
- ✓ Easily replaces existing 9 volt alkaline battery.
- ✓ Can be recharged without removal from Chair-Check II Control Unit.



Product No. 03-135-00200

Ohmmeter Connector

- ✓ Provides additional means of testing battery charge.
- ✓ Easy connection from Control Unit to ohmmeter.



Product No. 100126-101

Battery Recharger for Nickel-Metal Hydride Battery

- ✓ Can be left attached to Chair-Check II for use as an AC adaptor in stationary situations and keeps the nickel-metal hydride battery fully charged.
- ✓ Recharges nickel-metal hydride battery in as little as 14 hours.
- ✓ U.L. listed in the U.S. and Canada.
- ✓ Input: 120V AC, 60Hz
- ✓ Output: 10.5 V DC, 10mA



Product No. 100102-101

Wristbands



Product No. 304

Stickers



Product No. 307

Magnets



Product No. 932

"W" Adaptor for Bed-Check® Control Unit

- ✓ Convert a single nurse call room outlet to a multi-port outlet.



Product No. 100132-101

Universal Nurse Call "Y" Adaptor for Bed-Check Control Unit

- ✓ Allows use of nurse call button and Bed-Check unit from the same plug.
- ✓ Universal 1/4 inch phone plug.
- ✓ Pull-out protection sets off an alarm if cables are pulled out of nurse call adaptor.



Product No. 03-043-00100

Contact Bed-Check for all other Nurse Call Adaptor needs.

Bed-Check Control Unit Wire Holsters (Patent Pending)

- ✓ Optional bed or wall mounting for Model II, Model V and Model Vr Control Units.
- ✓ Optional chair or rail mounting for Chair-Check II Control Unit.
- ✓ Increase service life, reduce repairs.
- ✓ Easy monitor access for caregiver.



Product No. 100034-101 (1" hanger - bed)
Product No. 100034-102 (2" hanger - bed)

Product No. 100035-101 (1" hanger - chair)
Product No. 100035-102 (2" hanger - chair)

**ACCEPT NO
SUBSTITUTE
FOR THE
BED-CHECK
SYSTEM**

Bed-Check® products are developed, manufactured and supported by qualified engineering, biomedical and a trained medical sales force throughout North America.

1-800-523-7956 • www.bedcheck.com

FLEXIBLE SOLUTIONS FOR HEALTHCARE CHALLENGES



ACCESSORIES