

## Seat Belt Alarm Model 0707-945



### Features

- ▲ Bed-Check control unit alarms when Velcro® closure is undone
- ▲ Adjustable strap offers a range of 18" to 48"
- ▲ Fits most wheelchairs
- ▲ Rugged construction
- ▲ Compatible with Chair-Check® II and Basic-Check® units



### A great way to reduce the risk of falls

The Bed-Check® Seat Belt Alarm is designed to be used on wheelchairs to alert you if someone tries to exit the wheelchair. When the buckle is undone, the Bed-Check control unit (sold separately) will alarm. The alarm message can be configured in multiple ways with all of the advanced features of Bed-Check control units, such as multiple alarm tones.

**STANLEY®**

**Healthcare Solutions**

800-523-7956 | Fax 402-475-4023 | [www.bedcheck.com](http://www.bedcheck.com)