

IT Custom-App Sensormat®

Models 73040, 74040 & 74041



Features

- ▲ **Restraint-free**
- ▲ **Maintain dignity and privacy**
- ▲ **Compatible with all Bed-Check control units**
- ▲ **IT Custom-App Model 73040**
 - 4 ft. cord
- ▲ **IT Custom-App Model 74040**
 - 5.5 ft. cord
- ▲ **IT Custom-App Model 74041**
 - 8 ft. cord
- ▲ **Measures 3.5" x 19"**
- ▲ **30-day limited warranty**

A Great Way to Reduce the Risk of Falls

All the great features of Bed-Check® monitors can be used with our IT Custom-App Sensormat pad, giving you the flexibility to select the monitoring method that works best for each of your residents. Choose our discreet Sensormat pads when you need to know when the resident is leaving a bed, chair or wheelchair.

The IT Custom-App Sensormat pad has an intermediate length (19") that makes it suitable for a variety of uses: gurneys, bariatric wheelchairs, side chairs, bolster side mattresses, and others. The Sensormat pad is connected to a Bed-Check control unit and the individual is situated comfortably on top of the Sensormat pad. When the individual moves from the Sensormat pad, the Bed-Check control unit alarm sounds either with a local alarm or at the nurse call or both.

Distributed by:

STANLEY®

800-523-7956 | Fax 402-475-4023 | www.bedcheck.com

Healthcare Solutions